



**€13.950**  
plus tax and shipping

## A CLASSIC FLOAT TANK MADE IN SWEDEN

This affordable sensory deprivation tank has been used in Scandinavia since 2004. It's as big as a car, made of fibreglass, and has all of the hardware and features you'll need to achieve a perfect float environment at home or at your float spa.

Floatingrest International AB  
Bergshyddebacken 5  
S-187 76 TÄBY, SWEDEN  
info@floatingrest.com  
+46 70 715 76 62



# AQUA ANIMA

Floatation is a proven therapy for stress relief, relaxation and pain reduction.

# FLOAT TANK SPECS

## DIMENSIONS

**2.76 METERS LONG**

**1.50 METERS WIDE**

**1.30 METERS HIGH**

## LOAD ON SITE

**FLOAT TANK: 150 KG**

**FILTRATION SYSTEM: 50 KG**

**EPSOM SALT: 475 KG**

**WATER: 750 LITRES, KG**

## SANITATION

### CARTRIDGE FILTER

**H2O2 (Hydrogen Peroxide)**

**UV (Ultraviolet Light)**

The Aqua Anima float tank uses a magnetic drive centrifugal pump, a filter system and a UV water purification system. You can choose between the standard filter system with four smaller filters or a single Big Blue brand filter.

The ultraviolet (UV) lamp is effective in killing almost all bacteria and the manual addition of hydrogen peroxide (H2O2) increases and prolongs the effect of the UV treatment. The entire system is extremely user friendly, and it's easy to change the filters and the UV lamp.

## ROOM REQUIREMENTS

We recommend a 4 x 4 meter or a 3 x 5 meter space to fit the float tank and shower. You may wish to build a custom closet area to hide and protect the filtration system or place it behind a partition. Feel free to contact us about your particular circumstances.

If you are installing above ground floor, please be sure to check with the building owner about the load restrictions. The load of the float tank filled with water and salt is about 1400 KG or 3000 LBS. If in doubt, you may need to hire a structural engineer to assess whether the weight of the float tanks can be supported by the floors and walls.

The electrical system of the tank is hard-wired into the room (meaning there's no power cable which goes to an outlet) and operates in the following parameters: 230 Volts, 2.300 Watts and 10 Amps.

Want to discuss the different ways of fitting a float tank at your location? Contact us for a free assessment.

# SALT AND SET UP

**Epsom salt** is what makes you float. The salts are not included in the price of purchase, but we offer a Starter Pack (€600 per unit) which can be delivered with your order.

We will also assist you in locating a regional distributor of high-quality Epsom salts for future supplies. The salt is changed every 9 to 12 months at commercial float centres and much less frequently if the tank is used in a home.

The **Starter Pack** includes the following:

- + 26 kg hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)
- + 500 kg high quality Epsom salt (
- + 200 pairs of earplugs
- + Floating cushion to support stiff necks
- + Dosage sticks for measuring H<sub>2</sub>O<sub>2</sub>

**Professional installation** is available at no charge for labor though charges will apply for one installer's travel and lodging. PA pro install typically takes 3 days per unit to complete, and includes hands-on training with the equipment.

An **installation guide** and **operations manual** are included with your order if you wish to install the tank yourself. Phone, Skype and email support is also available for customers who choose to install the equipment by themselves.



# ABOUT

## MAINTENANCE

It is mandatory that you have a shower facility on the premises, and we recommend placing a shower as close as possible to the floatation units. This is for client comfort and for ease of maintenance because when the salt water drips off a client's body on the way to the shower, salt crystals form in the area of the drip. Much of the maintenance involved in operating a float center is in the manual cleaning of the float tank's interior and exterior, the shower area and filtration unit. Minimizing clean up will save you time, energy and frustration.

## BUSINESS MODEL

Two tanks should be invested in at the same time, or at least planned for at your location in case you can only purchase one now. The reason for this is a social one: it is common for two people to come together to float during the same time period. Couples, friends, a regular customer introducing a new person to the treatment ... these are typical examples. When only one device is available, people are deterred from coming together due to the wait time for their turn.

## WHAT IS FLOATING? THE BENEFITS

Floating Rest (Restricted Environmental Stimulation Technique) is a relaxation technique known in more than 30 countries throughout the world, and which has over 50 years of scientific research behind it.

Quite simply, you float upon water which is set to the body's temperature, and in an environment with subdued light, sound, movement and distraction. The water contains 30% Epsom salt (magnesium sulphate), a naturally occurring mineral compound known for its therapeutic qualities for skin and muscle.

A sensation of weightlessness is felt as your skin, the water and the air blend together in a sensory continuum. The point of floating is to relax, unwind and receive relief from pain. Floating works wonders on whiplash, rheumatoid arthritis, fibromyalgia and other conditions where people suffer from chronic stiffness, soreness and sensitivity. These are just some of the reasons for floating but it is especially attractive to patients because no pills, injections or operations are involved. Sensory isolation is also used for super learning, spiritual development and mental training for athletes. It is worth first giving the brain the chance of repairing damages caused for different reasons. Remember the brain uses 85% of its capacity just keeping the body balanced, operational and aware of its environment. It is easy to understand the power of this resource when it is freed to troubleshoot in the body and mind. That's when the "miracles" occur and when people fall in love with floating.

**NEED MORE INFO?**

**+46 70 715 76 62**

**OR REQUEST A  
CALLBACK**

[WWW.FLOATINGREST.COM](http://WWW.FLOATINGREST.COM)

[info@floatingrest.com](mailto:info@floatingrest.com)

*Floatingrest International AB  
Bergshyddebacken 5  
S-187 76 TÄBY, SWEDEN*